



Main technical parameters

Implementation standards

• Temperature fluctuation: ≤±0.5 °C

GB/T5170.2-2008 Temperature test equipment

Temperature uniformity: ≤2.0 °C

• Temperature deviation: ≤±2.0 °C

GB/T2423.1-2008(IEC68-2-1) testing A, Low temperature test method

Ambient temperature: +5~+35°C

Power(V): AC 380±10%V 50HZ±0.5HZ

• GB/T2423.2-2008(IEC68-2-2) testing B, High temperature test method

Ambient temperature: +5~+35°C


Equipment noise: ≤69dB

• GJB150.3A-2009(MIL-STD-810F-2000) High Temperature test

(testing from one meter in f/ront of the door)

- Standard configuration: Electrothermal film glass observation 2 pcs;
- GJB150.9A-2009(MIL-STD-810F-2000) thermal humidity test (C)

Standard configuration: Electrothermal film glass observation 2 pcs; Cable hole (Φ 100) 1 PCS; sample shelf 2 kits; Lighting 1 pcs; Sample power control terminal 1 (C), only C type equipment equipment with this.



Model		SH-222	SH-242	SH-262	SH-642	SH-662	SH-242-5	
Syste	em		Balanced Temperature & Humidity Control system (BTHC system)					
Temp. performance ^{*1}	Temp. range	- 20 to +150°C (- 4 to +302°F)	- 40 to + 150°C (- 40 to + 302°F)	- 60 to + 150°C (- 76 to + 302°F)	- 40 to + 150°C (- 40 to + 302°F)	- 60 to + 150°C (- 76 to + 302°F)	- 40 to + 150°C (- 40 to + 302°F)	
	Temp. fluctuation	±0.3°C (-20to +100°C) ±0.5°C (+100.1 to +150°C)	±0.5℃	±0.3°C (-60to +100°C) ±0.5°C (+100.1 to +150°C)	±0.3°C (-40to +100°C) ±0.5°C (+100.1 to +150°C)	±0.3°C (-60to +100°C) ±0.5°C (+100.1 to +150°C)	±0.3°C (-40to+100°C) ±0.5°C (+100.1 to +150°C)	
	Temp. gradient / Temp. variation in space	2.5°C (-20to +100°C) 4.0°C (+100.1 to +150°C)	2.5°C (-40to +100°C) 4.0°C (+100.1 to +150°C)	2.5°C (-60to +100°C) 4.0°C (+100.1 to +150°C)	2.5°C (-40to +100°C) 4.0°C (+100.1 to +150°C)	2.5°C (-60to +100°C) 4.0°C (+100.1 to +150°C)	2.5°C (-40 to +100°C) 4.0°C (+100.1 to +150°C)	
	Temp. Heat up rate		3.2℃ /min.		2.9℃	/min.	5.0℃ /min.	
	rate of Pull down rate		2.1℃ /min.		1.7℃	/min.	5.0°C /min.	
	Temp. extreme achievement time Heat up time	From - 20 to + 150℃ within 55 min.	From - 40 to +150°C within 60 min.	From - 60 to +150℃ within 70 min.	From - 40 to +150°C within 70 min.	From - 60 to +150℃ within 80 min.	From - 40 to + 150℃ within 40 min.	
	Temp. extreme achievement time Pull down time	From + 20 to - 20°C within 20 min.	From + 20 to - 40°C within 50 min.	From + 20 to - 60°C within 70 min.	From + 20 to - 40°C within 60 min.	From + 20 to - 60°C within 90 min.	From + 20 to - 40°C within 20 min.	
	Lowest attainable tem	p 20℃	- 40℃	- 60℃	- 40℃	- 60℃	- 40℃	
Humid. performance ^{*1}	Humid. range	umid. range 30 to 95% rh (Refer to diagram on page 12)						
	Humid. fluctuation ±3.0% rh							
Construction	Heater Nichrome strip wire heater							
	Humidifier							
	System	Mechanical single-stage refrigeration system Mechanical cascade refrigeration system						
	Cooler	Plate fin cooler						
	Refrigerator Refrigerator capaci			(, Air-cooled condenser, Expansion mechanism: Capillary tube system [Unit 1: 400W ×1, Unit 2: 400W ×1]				
	∝ Refrigerant	R404A		R23,R404A				
Capa	acity		22.5 L		64	4 L	22.5L	
	mber total load stance	20 kg						
Insid	le dimensions (inch) *2		W300×H300×D250 (W11.81×H11.81×D9.84		W400×H4 (W15.75×H15	00×D400 .75×D15.75)	W300×H300×D250 (W11.81×H11.81×D9.84)	
	Outside dimensions W440×H6 mm (inch) *2 (W17.32×H2)		90×D695 (18×D27.36)	W440×H690×D785 (W17.32×H27.18×D30.91)		/30×D920 3.74×D36.22)	W440×H690×D785 (W17.32×H27.18×D30.91)	
Weig	/eight 83 kg (78 for 100		r 100V type)	105 kg 130 kg		106 kg		
nts	Allowable ambient conditions		+ 5 to + 35°C (+ 41 to + 95°F)					
	100V AC 1 φ50/60Hz		5 A	18.0 A	21.0 A	.0 A	21.0 A	
115V AC 1 φ60Hz 200V AC 1 φ50/60Hz *4		14.0 A						
				14.0 A	14.	.5 A	15.5 A	
	220V AC 1 φ50/60Hz *5	10.	0 A	13.5 A	14.	0 A	15.0 A	
230V AC 1 ф50Hz ^{*5}		9.5	9.5 A		14.	0 A	15.0 A	
Noise level *6		Between 42	Between 42 and 52 dB		Between 4	8 and 59 dB	Between 42 and 56 dB	
Exhaustheatquantity		3500 kJ/h(836 kcal / h)	4000 kJ/h (955 kcal/h)	5040 kJ/h (1204 kcal/h)	5700 kJ/h (1361 kcal/h)	

CLIMATE STSAR series has advanced features in terms of quality and reliability

Customer first

- 1. If you have ever used environmental test equipment, you will soon feel the unique design and ease of use of the device CLIMATE STSAR.
- 2. First of all, you can feel the equipment is easy to use, low maintenance rate and high reliability
- **3.** Then, You can choose different the testing volume, temperature range and special parts to meet your special requirements


Products Features

The CLIMATE STAR series products have excellent design and high quality standard features.

- 1. Large viewing angle and full heating observation window
- 2. High stability full color touch screen
- 3. Pin holes on both sides
- 4. Sample holder capable of conveniently adjusting height
- 5. Triple independent over temperature protection
- 6. Safety sample terminal
- 7. High quality casters for the equipment easily level shift
- 8. Ultra quiet
- 9. Disassemble operation panel for easy maintenance

Scope of application

- 1. This model is a simulation products in climate field, it's combination of temperature conditions (high and low temperature operation & storage, temperature cycle, high temperature, low temperature, condensation test) testing the product whether it has any changes in the ability and the characteristic.
- 2. Must meet the requirements of the international standard (IEC, JIS, GB, MIL---) to achieve the consistency of the international



SM C 1000 C


Humidifying laver

D

Indicate temperature range

Programmable

Internal volume of the equipment
Mild (Energy saving)
5Brand code/name(sanwood)



Structure characteristics

- 1. Shell: Spray galvanized color steel plate, the surface electrostatic spray processing.
- 2. Liner: stainless steel SUS 304.
- 3. Thermal insulation layer: Polyurethane foam board and glass fiber.
- 4. Seal: Toshiba high purity silicon rubber raw materials, effectively prevent aging.
- **5.**Heater: Ni Cr alloy electric heater.
- 6. Humidifier: Outer tube: SUS316 stainless steel seamless pipe Internal heating wire: Ni Cr alloy wire.
- 7. Sample holder: 40kg/ layer * 2 layer (standard configuration)
 80kg/ layer ; 120kg/ layer
 Total bearing ≤ 240 kg (optional).

Pin hole



Φ100mm (standard equipment) Φ50mm Φ80mm Φ160mm (optional equipment)

Operation sample hole on the glass (No display on the picture)



Inner glass door (optional)



Refrigeration design

- 1. Modular production, reliable quality, convenient maintenance.
- 2. Silver brazing welding vibration pipe with a silver content of 45% to prevent the welding leak effectively.
- 3. Adequate space position, easy to operate.
- 4. Welding through nitrogen, ensure the inner pipe not nitriding.
- 5. Take a variety of techniques to decouple shock.
- 6. Take a variety of techniques to anti-corrosive.



Compressor







France Taikang compressor (Original import)

Pressure relay



America EMERSON or Denmark DANFOSS

Evaporator

Custom efficient fin type heat exchanger

Solenoid valve



Italy CASTEL

Refrigerant

R404A R23(-70) Ozone depletion index was 0

Denmark DANFOSS brand



- 1.condenser
- 2.evaporator condenser(-70)
- 3. Evaporation pressure regulating valve
- 4. Thermal expansion valve
- 5.Dry filter
- 6.Condensation pressure regulating valve (water-cold)

Controller



- 1. 5.7" 640*480 lattice. TFT LCD displayer
- 2. 1200 programs, program can cycle
- 3. RS 485 interface, with remote communication function.
- 4. SD card storage test data, about7500 days (Sampling period: 5min)
- 5. operating language: Chinese or English

Recorder(option)



- 1. Large screen LED display
- 2. High reliability of industrial r ecords requirements

The sample power control terminal



1. When the equipment safety protection device works, the power supply of the electrified sample is controlled through the connecting terminal.

Safety protection device

1.Compressor



- 1.1 Compressor overpressure
- 1.2 Compressor motor overheating
- 1.3 Compressor motor over-current
- 1.4 Condenser fan overheating (air-cold)
- 1.5 Cooling circulating water pressure shortag (water-cold).

2. Test samples of protection



- 2.1 Adjustable overtemperature protection.
- 2.2 Air conditioning channel over temperature limit.
- 2.3 controller set overtemperature shut down alarm.
- 2.4 sample terminal protection.

3. Electric control

- 3.1 The fan motor overheating.
- 3.2 Total power phase sequence and lack of phase protection.
- 3.3 Leakage protection.
- 3.4 Load short circuit protection.

The Experience you Rely on...

Sanwood Environmental Chambers was established in 1995, which integrated Taiwan and Japan technologies. We have been focus on the most secure and reliable climatic test chamber technology since established. And has become a private science and technology enterprises in Dongguan, Guangdong Province, which passed the ISO9001:2008 quality system certification.

Our products upgrade constantly and our customers come portable batteries, power batteries, battery, lithium batteries, lead-acid, new energy vehicles, electric bicycles, electric tools, electric systems, solar, military, universities research and other technology industries fields.

Having experienced nearly 20 years efforts, we have successfully developed a series of products:

power batteries, battery, lithium batteries, lead-acid, new energy vehicles, electric bicycles, electric tools, electric systems, solar, military, universities research and other technology and stries applies applies a the technology and stries applies are the technology and stries are the technology and the technology are the technology are the technology and the technology are the technology and the technology are the technology are the technology and the technology are the technology are the technology are the technology and the technology are the technol

- an explosion-proof type temperature test box
- walk-in temperature and humidity chamber
- weather resistance test chamber
- battery thermal abuse test box
- explosion-proof type h ot box
- Tem perature&hum idity&Vibration integrated test chamber
- du st test box
- vib ration table
- rain test chamber
- oz one test box
- xenon lamp test chambe
- high temperature oven
- seaw ater immersion box

All of products meet GB31241、IE62133、QCT/743、UN38.3、UL2054 Standard. And we have had a good cooperation with ATL, Sony, Sunwoda, Desay, Samsung, BYD, Toyota, Yutong Bus, Nissan, Guangdong Province entry-exit, Tsinghua University, Henan University, Chinese Academy of Sciences, Central South University Successively.

Enterprise vision:

Sanwood Technology has established a large production base in Dongguan after many years efforts. The plant area reached more than 12000 square meters. The foreign trade branch and foreign service agencies were established in 2010. And branches successively established in Taiwan, Suzhou, Hunan, Hubei, Beijing, Henan. Excellent products and good after-sales service make us won the recognition and trust of customers. Products are exported to more than 30 countries, such as Russia, Singapore, the United States, Turkey, Denmark, Vietnam, India, Malaysia, Kazakhstan, Austria, Canada, etc. In the age with fierce competitions, Sanwood thrived little by little and aims to become the leading brand in the safety and reliability environmental test













Focusing on the innovation of environmental reliability test

Sanwood(HK)Industrial Corporation.,Limited
Guang dong Sanwood Technology Corporation.,limited

ADD:Changping Science&Technology park,Changping, Dongguan city,Guangdong province,China

TEL:+86+769-81182799

FAX:+86+769-82987199

E-mail:info@climatic-chambers.com.tw

www.climatic-chambers.com.tw

www.environmental-chambers.ru